





## PCR Tube Racks 96-Well

### Perfect for pre- or post-PCR sampling

- Made of polypropylene
- Accommodates individual tubes or strips
- Removable hinged lid lies flat
- Alphanumeric identification for quick location and retrieval of tubes



Pack of 5

Item No.	Color	Item No.	Color
HS23461A	Assorted* 	HS23461C	Blue 
HS23461B	Natural 	HS23461D	Green 

\*Assorted colors (Blue, Green, Pink, Yellow, Orange)

Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
0.2	96	5.2	5 x 3.9 x 1.2	12.7 x 9.9 x 3.0



**CORPORATE OFFICE**  
**GENAXY SCIENTIFIC PVT. LTD**

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058  
Ph: +91-11-4761 9999(100 Lines) Web: [www.genaxy.com](http://www.genaxy.com) E-Mail: [technical@genaxy.com](mailto:technical@genaxy.com)