

PCR Tube Racks Reversible



Removable hinged lid fits both sides

- Made of polypropylene
- PCR side of the rack holds 8 or 12 x 0.2 mL tube strips or individual tubes
- Reverse side holds 0.5 mL tubes
- Both sides of the rack have 12 wells for 1.5 mL tubes

Pack of 5

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in (cm)
HS2344A		0.2	168	5.5	8.1 x 5.2 x 2.1 (20.6 x 13.1 x 5.3)
		0.5	40	7.6	
		1.5/2.0	12	11.4	

**Assorted colors (Blue, Green, Pink, Yellow, Orange)*



Engineering, OEM & Industrial Design

- We have in-house product engineering and R&D teams that are dedicated to delivering quality, innovative, efficient, and compact solutions that customers can rely on.
- We have a range of capabilities from private labeling to product development, design, packaging and manufacturing.



CORPORATE OFFICE
GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com