


Low-Temp PCR Rack



Sliding Lid
rotates 180° to become a convenient stand the rack sits on

US Patent Pending, Community Design No: 004101756-0001-0002

Pack of 5

Item No.	Color	Item No.	Color
120538	Assorted*	120540	Green
120539	Blue	120541	Purple

*Assorted colors (Purple, Blue, Green, Orange, Yellow)



Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
0.2	96	5.2	4.9 x 3.6 x 1.5	12.5 x 9.2 x 3.7

Low-Temp PCR Rack: 96-Well

Slide, set, work & store PCR samples

Multi-purpose rack converts from compact workstation to lidded storage solution accommodating individual and strip 0.2 mL tubes

Versatile & Robust Solution

- Polycarbonate lid serves as cover and elevated rack base - just slide open and pivot lid under rack for compact workstation
- Lid can be used in right or left-handed orientation and is detachable
- Stackable saving valuable bench-top and freezer space
- Multiple colors offered for coding options

CORPORATE OFFICE

GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058

Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com