


Reversible Tube Racks 96-Well




Tabs on each end of rack snap lid into place

Raised Rings around the 1.5 to 2.0 mL wells make tubes easier to remove



Molded Grid Lines

0.5 mL wells separated by molded-in grid lines, making them easier to see



Reversible Tube Racks 96-Well

Easy-to-grasp tabs for sample handling and storage

96 wells for 0.5 mL tubes on one side and 96 wells for 1.5/2.0 mL tubes on the reverse side

Optimum Function & Identification

- Wells arranged in standard 8 x 12 array
- Molded in-grid lines separate 0.5 mL wells making them easier to see
- Raised rings around the 1.5 mL/2.0 mL wells make tubes easier to remove
- Alphanumeric identification for quick location and retrieval of tubes

Secure lid

- Translucent lid snaps onto rack tabs for secure, visible storage
- Spacious enough to fit most screw-top micro-centrifuge tubes
- Stackable, saving valuable space



Pack of 5

Item No.	Color	Item No.	Color
HS2345A	Assorted*	HS2345D	Green
HS2345B	Natural	HS2345E	Pink
HS2345C	Blue	HS2345F	Yellow

*Assorted colors (Blue, Green, Yellow, Orange, Pink)

Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
1.5/2.0	96	11	9.6 x 4.7 x 2	24.4 x 12 x 5.1
0.5		7.6		

CORPORATE OFFICE
GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com