


Rota-Rack® Duo Tube Rack

Easy Identification

Alphanumeric well identification



Rota-Rack® Duo Tube Rack

Dilution and concentration studies made easy

Transfer solutions between 0.2 mL PCR tubes and 0.5 and 1.5 mL microcentrifuge tubes

- Purple section holds 48 x 0.5 mL and 48 x 1.5/2.0 mL microcentrifuge tubes on the reverse side
- Blue section holds 96 x 0.2 mL tubes and 48 x 1.5/2.0 mL tubes on the reverse side



Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS29041A	Blue/ Purple	0.2	96	5	11.6 x 4.6 x 1.3	29.6 x 11.8 x 3.4
		0.5	48	8		
		1.5/2.0	96	11		

CORPORATE OFFICE
GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com