

Microtube Racks 80-Well

Popular fraction-collector-plate style rack

Compact 5x16 array makes filling and transfer easy












- Frosted labeling area on bottom of rack
- Embossed index numbers for quick sample ID and retrieval
- Made of polypropylene

Tube Racks & Storage



Pack of 5



Item No.	Color	Part No.	Color
HS29025A	Natural	HS29025E	Yellow 
HS29025B	Blue 	HS29025F	Orange 
HS29025C	Green 	HS29025G	Assorted*     
HS29025D	Pink 	HS29025J	Red 

*Assorted colors (Blue, Green, Yellow, Orange, Pink)

Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
1.5/2.0	80	11	8.9 x 2.6 x 1	22.5 x 6.7 x 2.8

CORPORATE OFFICE
GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com