

Microtube Racks

Durable, Simple to Use

- Made from reinforced nylon offering superior chemical and heat-resistance over most other plastic based racks
- Handles provide easy and safe transportation
- Stable design can be used in water baths
- Alphanumeric identification for quick location and retrieval of tubes



Item No.	Color	Tubes mL	Wells	Well Dia. mm	Packaging	L x W x H in	L x W x H cm
120241	Blue ■	1.5/2.0	24	11.1	8	5.6 x 3.5 x 2	14.3 x 9 x 5.2
120243			96		4	8.3 x 5.2 x 2.1	21 x 13.1 x 5.4

CORPORATE OFFICE
GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com