


## 81-Well Microtube Storage Boxes

### Microtube Storage Boxes

- Hold 1.5 and 2.0 mL microcentrifuge tubes and cryogenic vials



### Alphanumeric Identification

for quick location and retrieval of tubes

### Stylish & Identifiable

- Made of translucent, vibrantly colored polypropylene allows for visibility inside box and color coding identification
- Tube capacity is indicated on the outside the box within a distinctive wave design
- Quick visual identification with uniform wells and alphanumeric location identification along with a frosted labeling area

### Robust

- Easy open and secure lid fitted with either friction fit or snap clasp designs
- Stackable, saving valuable space in the lab
- Black, light blocking version available for light sensitive samples



Pack of 5

Item No.	Color	Item No.	Color
120034	Natural	120038	Yellow
120035	Blue	120039	Orange
120036	Green	120040	Assorted*
120037	Pink	120041	Black

\*Assorted colors (Blue, Green, Yellow, Orange, Pink)

Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
1.5/2.0	81	12.3	5.1 x 5.1 x 1.9	13 x 13 x 4.7

**CORPORATE OFFICE**  
**GENAXY SCIENTIFIC PVT. LTD**

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058  
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com