



## 100-Well Microtube Storage Boxes

### Microtube Storage Boxes

- Hold 1.5 and 2.0 mL microcentrifuge tubes and cryogenic vials



### Alphanumeric Identification


for quick location and retrieval of tubes

### Stylish & Identifiable








- Made of translucent, vibrantly colored polypropylene allows for visibility inside box and color coding identification
- Tube capacity is indicated on the outside the box within a distinctive wave design
- Quick visual identification with uniform wells and alphanumeric location identification along with a frosted labeling area

### Robust

- Easy open and secure lid fitted with either friction fit or snap clasp designs
- Stackable, saving valuable space in the lab
- Black, light blocking version available for light sensitive samples



Pack of 5

Item No.	Color		Item No.	Color	
120042	Natural		120202	Green	
120043	Black		120203	Yellow	
120044	Assorted*		120204	Orange	
120201	Blue				

\*Assorted colors (Purple, Blue, Green, Yellow, Orange, Pink)

Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
1.5/2.0	100	12.3	5.6 x 6 x 2.2	14.1 x 15.1 x 5.7

### CORPORATE OFFICE GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058  
Ph: +91-11-4761 9999(100 Lines) Web: [www.genaxy.com](http://www.genaxy.com) E-Mail: [technical@genaxy.com](mailto:technical@genaxy.com)