


Rota-Rack® Tube Racks

Get organized with this interlocked rack

Each module is interlocked and adjustable with just a twist allowing you to work with multiple sized tubes in the same rack

- Suitable for variety of culture, centrifuge, micro-centrifuge, round, flat or conical bottomed tubes
- Two versions to accommodate small and large tube sizes




Rotate and Customize
to suit your application needs

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS29040A	Yellow	50	2	29	8.8 x 4.3 x 4.3	22.4 x 11 x 11
	Green	10-15	8	20		
	Blue	15-20	6	16		
	Purple	10-13	10	13.5		
HS29040B	Yellow	15	6	16.4	8.2 x 3.9 x 3.9	21 x 10 x 10
	Green	1.5/2.0	9	11.2		
	Blue	0.5	12	8		
	Purple	0.2	32	6.5		

Version available for PCR and Microtubes:

Rota-Rack® Duo Tube Rack



CORPORATE OFFICE
GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com