

Cube Tube Rack

Interlocking multi-tube racks

- Made of polypropylene
- Allows you to work with different sized tubes all in one unit



Pack of 5

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS29050A		50	4	30	4.3 x 4.3 x 4.5	10.8 x 10.8 x 11.3
		15	10	17		
		3-7	12	12.5		
		1.5/2.0	16	11.5		

*Assorted colors (Blue, Green, Pink, Yellow, Orange)

CORPORATE OFFICE
GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com