

4-Way Tube Rack



4-Way Tube Rack

Link together to meet your application needs

Versatile

- Use as standalone unit or infinite units can be linked to create configurations to fit number of wells and tube sizes needed

Robust Functional Design

- Molded alphanumeric identification for quick location and retrieval of tubes
- Made from polypropylene that is chemical resistant and easy to clean
- Racks can be used in water baths and placed in freezer



Pack of 5

Item No.	Color	Item No.	Color
HS29022A	Natural	HS29022E	Yellow
HS29022B	Blue	HS29022F	Orange
HS29022C	Green	HS29022G	Assorted*
HS29022D	Pink	*Assorted Colors (Blue, Green, Pink, Yellow, Orange)	

Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
50	4	30	6.9 x 3.7 x 2.1	17.4 x 9.5 x 5.2
15	12	17.5		
0.5	32	8.2		
1.5/2.0	32	11.6		

CORPORATE OFFICE GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com