

Transformer Cube Tube Rack

Low temp storage of large samples,
15 mL & 50 mL tubes

- Quick change over with 2 snap-in inserts that hold 15 and 50 mL tubes
- Rack keyed and numeric wells for tube identification
- Clear Lid included for viewable, stackable storage
- Open space available under insert for extra storage



Pack of 5

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H Dia.	L x W x H cm
HS23271A	Purple	15	16	18	4.9 x 4.9 x 5	12.5 x 12.5 x 12.9
		50	9	31		

CORPORATE OFFICE
GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com