

## OneRack® Multi

### Exclusive design, optimal functionality

The most flexible tube rack designed for use with a multitude of tube diameters and heights

Holds virtually any brand tube 10-30 mm of varying heights


### Multiple Configuration for Workflow Efficiency

- Interchangeable tube rack inserts hold multiple tube sizes firmly in place when alternating between sample vessels
- Inserts can be pivoted in either an upper or lower position to securely accommodate the size and height of tube being used


### Multi-Functional Design

- Open design allows for better filling visibility of tube graduations when filling
- Rack wells have a unique teardrop design that enables better retention of tubes over traditional designs
- Rotate tube rack insert to an angle position to facilitate more comfortable pipetting

Unique Well Design  
Tear drop shape allows for better retention of tubes




Lower position is ideal height when working with shorter microcentrifuge tubes



Alphanumeric grid for quick sample identification and retrieval

Angled position for more comfortable pipetting



Ergonomic handles for ease of transporting the rack



Item No.	Color	L x W x H in	L x W x H cm	# Inserts Included	Well Dia./Insert	Wells/Insert
120817	Blue/Natural	12.2 x 4.2 x 3	30.9 x 10.7 x 7.7	6 (2 ea/ Diameter)	13 mm 17 mm 30 mm	20 20 6
120818	Purple/Natural	12.2 x 4.2 x 3	30.9 x 10.7 x 7.7	6 (2 ea/ Diameter)	13 mm 17 mm 30 mm	20 20 6

### OneRack® Multi

US Patent No. 9776188

Item No.	Color	L x W x H in	L x W x H cm	# Inserts Included	Well Dia./Insert	Wells/Insert
120817	Blue/Natural	12.2 x 4.2 x 3	30.9 x 10.7 x 7.7	6 (2 ea/ Diameter)	13 mm 17 mm 30 mm	20 20 6
120818	Purple/Natural	12.2 x 4.2 x 3	30.9 x 10.7 x 7.7	6 (2 ea/ Diameter)	13 mm 17 mm 30 mm	20 20 6

### CORPORATE OFFICE GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058  
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com