

## One-Well Connecting Tube Racks

Form a line or use as a single, stand-alone polypropylene rack

- Open Design enables easy identification of samples
- Molded circular base ring keeps conical or round-bottom tubes in place

Pack of 5

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS23052		50 *	1	29	3.2 x 2.8 x 2.6	8.3 x 7 x 6.7

\*Assorted colors (Blue, Green, Yellow, Orange, Pink)



## One-Well Friction-Fit Tube Rack

- Made from polystyrene
- Open Design enables easy identification of samples
- Holds conical and round bottom tubes
- Slides smoothly across the benchtop

## One-Well Wire Tube Rack

- HDPE-coated steel wire
- Holds conical and round bottom tubes
- Open Design enables easy identification of samples
- Base ring keeps tubes in place



### One-Well Friction-Fit Tube Rack

Pack of 5

Item No.	Color	Tubes mL	Wells	Well Dia. mm	W x H in	W x H cm
HS23053A	Blue	50	1	29	3.6 X 1.2	9.2 x 3

### One-Well Wire Tube Rack

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS23051	Yellow	50	1	29	3.5 x 2.6 x 2.6	9 x 6.7 x 6.7



**CORPORATE OFFICE**  
**GENAXY SCIENTIFIC PVT. LTD**

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058  
Ph: +91-11-4761 9999(100 Lines) Web: [www.genaxy.com](http://www.genaxy.com) E-Mail: [technical@genaxy.com](mailto:technical@genaxy.com)