

## Snap-Together Conical Tube Rack

Space Saver - shipped flat and easy-to-assemble

Conveniently holds 15 and 50 mL tubes simultaneously

- ABS racks snap together securely and are designed for stable use in water baths
- Alphanumeric identification for quick location and retrieval



Item No.	Color	Tubes <i>mL</i>	Wells	Well Dia. <i>mm</i>	L x W x H <i>in (cm)</i>
HS24306	Orange <span style="color: orange;">■</span>	15	30	17	8.1 x 6.7 x 2.4
		50	20	28	(17.1 x 17 x 6.1)

**CORPORATE OFFICE**  
**GENAXY SCIENTIFIC PVT. LTD**

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058  
Ph: +91-11-4761 9999(100 Lines) Web: [www.genaxy.com](http://www.genaxy.com) E-Mail: [technical@genaxy.com](mailto:technical@genaxy.com)