



OneRack® Tube Rack

Tough and highly resistant to acids, bases, solvents and heat

Durable racks are molded in a single, continuous piece so no assembly required and no detachable pieces

Intuitive Design

- Unique tear drop shaped openings offer better retention of tubes
- Molded in coordinates for tube orientation
- Flat labeling area for easier identification
- Stable design for use with water baths

Available in two materials:

Fiberglass-Reinforced Polyoxymethylene (POM)

- Available in half and full rack sizes for 5 different tube diameters
- Offers the highest chemical and heat resistance (-40 to +121°C)

Polypropylene

- Available in full rack sizes for 2 different tube diameters
- Offers good chemical and heat resistance (0 to 110°C)

OneRack Full-Size POM



Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in (cm)
HS27511A	White				
HS27511B	Blue				7.9 x 4 x 2.3
HS27511C	Green	5-10	72	13	(20.2 x 10.4 x 5.9)
HS27511D	Magenta				
120205	Yellow				
HS27512A	White				
HS27512B	Blue				9.8 x 5 x 2.7
HS27512C	Green	5-10	72	16	(25 x 12.7 x 7)
HS27512D	Magenta				
120206	Yellow				
HS27513A	White				
HS27513B	Blue				9.9 x 3.9 x 3.2
HS27513C	Green	10-18	40	20	(25.2 x 10 x 8.3)
HS27513D	Magenta				
HS27514A	White				
HS27514B	Blue				11.8 x 4.7 x 3.6
HS27514C	Green	10-18	40	25	(30 x 12 x 9.2)
HS27514D	Magenta				
HS27515A	White				
HS27515B	Blue				11.1 x 4.3 x 3.3
HS27515C	Green	15-50	24	30	(28.2 x 11 x 8.5)
HS27515D	Magenta				
120207	Yellow				

OneRack Half-Size POM



Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in (cm)
HS27501A	White				
HS27501B	Blue				
HS27501C	Green	5-10	36	13	4.1 x 4.1 x 2.3 (10.4 x 10.4 x 5.9)
HS27501D	Magenta				
120398	Yellow				
HS27502A	White				
HS27502B	Blue				
HS27502C	Green	5-10	36	16	5 x 5 x 2.7 (12.7 x 12.7 x 7)
HS27502D	Magenta				
120399	Yellow				
HS27503A	White				
HS27503B	Blue				
HS27503C	Green	10-18	20	20	4 x 5 x 3.2 (10 x 12.7 x 8.3)
HS27503D	Magenta				
HS27504A	White				
HS27504B	Blue				
HS27504C	Green	10-18	16	25	4.7 x 4.8 x 3.6 (12 x 12.2 x 9.2)
HS27504D	Magenta				
HS27505A	White				
HS27505B	Blue				
HS27505C	Green	15-50	9	30	4.3 x 4.3 x 3.3 (11 x 11 x 8.5)
HS27505D	Magenta				

OneRack Full-Size Polypropylene

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in (cm)
HS27551A	White				
HS27551B	Blue				
HS27551C	Green	5-10	72	13	7.9 x 4 x 5.9 (20.2 x 10.4 x 5.9)
HS27551D	Magenta				
HS27552A	White				
HS27552B	Blue				
HS27552C	Green	5-10	72	16	9.8 x 5 x 2.7 (25 x 12.7 x 7)
HS27552D	Magenta				

CORPORATE OFFICE

GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com