

## Universal Tube Rack



### Adjustable

Twist the knob to adjust the well sizes for your tubes

### Ideal for holding a wide range of tubes

- Made of polypropylene
- Holds up to 12 similar size tubes
- Alphanumeric identification for quick location and retrieval of tubes
- Side handles for easy transport

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
<b>HS24330A</b>	White/Purple	□ ■ 15–50	12	16–30	7.6 x 5.5 x 2	19.4 x 14 x 5.1

**CORPORATE OFFICE**  
**GENAXY SCIENTIFIC PVT. LTD**

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058  
Ph: +91-11-4761 9999(100 Lines) Web: [www.genaxy.com](http://www.genaxy.com) E-Mail: [technical@genaxy.com](mailto:technical@genaxy.com)