


Fold and Snap Tube Rack

Economical durable rack

- Polypropylene racks ship flat in one piece
- Easy-to-assemble, just fold and snap together
- Stable design can be used in water baths
- Open design allows samples to be incubated at the same temperature
- Color code for easy identification



Easy to assemble
Just fold and snap together!



Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS243073B	Blue	5-10	90	13	9.6 x 4.1 x 2.4	24.5 x 10.5 x 6.3
HS243073W	White					
HS243073R	Red					
HS243073Y	Yellow	5-15	60	16	9.6 x 4.2 x 2.8	24.5 x 10.8 x 7.2
HS243077B	Blue					
HS243077W	White					
HS243077R	Red	5-30	40	20	9.8 x 4.2 x 2.8	24.9 x 10.8 x 7.2
HS243077Y	Yellow					
HS243071B	Blue					
HS243071W	White					
HS243071R	Red					
HS243071Y	Yellow					

CORPORATE OFFICE
GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com