


Mega Rack®

High-capacity storage for your tubes

Fully assembled, polypropylene mega racks are great for high throughput labs requiring large sample storage

- Sturdy, material for easy and secure transport of blood collection and other tubes
- Alphanumeric identification for quick location and retrieval of tubes
- Place different test groupings into one rack eliminating the need for multiple racks



Mega Racks® 10-13 mm - Holds 5-7 mL sample tubes

Item No.	Color	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS24322A	Blue ■	216	13	13.3 x 9.3 x 3.6	33.7 x 23.5 x 9.1
HS24322E	White 				
120112	White 	432	13	26.5 x 9.3 x 3.5	67.3 x 23.5 x 9.0
120113	Blue ■				

Mega Racks® 13-16 mm - Holds 10 mL sample tubes

Item No.	Color	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS24323A	Blue ■	120	16	13.3 x 10.5 x 3.6	33.7 x 26.7 x 9.1
HS24323E	White 				
120114	White 	240	16	26.5 x 10.5 x 3.5	67.3 x 26.7 x 9.0
120115	Blue ■				

Mega Racks® 18-20 mm - Holds 15 mL conical/swab tubes and 20 mm round bottom tubes

Item No.	Color	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120834	White 	120	20	14.9 x 11.9 x 3.6	37.7 x 30.1 x 9.1
120835	Blue ■				

CORPORATE OFFICE
GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com