



## Peg Tube Racks

Strong, tough, versatile racks for bench use or tube storage

Fiberglass-reinforced polypropylene racks are ideal for bench work or tube storage

- Spacious peg design allows easy viewing of tube contents
- Pegs remain straight even after repeated autoclaving
- Embossed alphanumeric identification for quick location and retrieval of tubes



Unique tear drop well design allows for better retention of tubes

Pack of 2

| Item No. | Color                                    | Tubes<br>mL | Wells | Well Dia.<br>mm | L x W x H<br>in | L x W x H<br>cm   |
|----------|--|-------------|-------|-----------------|-----------------|-------------------|
| HS24311A | White <input type="checkbox"/>           | 10          | 50    | 14-17           | 7.4 x 4.1 x 2.7 | 18.8 x 10.4 x 6.9 |
| HS24311B | Blue <input checked="" type="checkbox"/> |             |       |                 |                 |                   |
| HS24312A | White <input type="checkbox"/>           | 7           | 96    | 10-13           | 3.2 x 5 x 2.7   | 8.2 x 12.6 x 6.9  |
| HS24312B | Blue <input checked="" type="checkbox"/> |             |       |                 |                 |                   |