


Floating Microtube Rack

For use on bench and in ice and water baths

- Made of polypropylene
- Attach legs and rack stands upright for use on bench top
- Detach legs and samples can be quickly removed by pushing down onto any surface



Pack of 4

Item No.	Color	Tubes mL	Wells	Well Dia. mm	Shape/Dia. in
HS2132A	White	1.5/2.0	16	10.9	Square 4
HS2134A			20		Round 3.8

CORPORATE OFFICE
GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com