

Cryogenic Vial Rack


Embossed Index Numbers and squared, upper left corner provide easy orientation



Ribbed Grips make handling with gloves easier

Safely store and transport cryogenic vials

- Made from polycarbonate
- Well tabs allow one-handed opening and closing of vials that have a slot feature
- Embossed index numbers and squared, upper left corner provide easy orientation
- Drainage holes in each well aid in cleaning
- Alphanumeric identification for quick location and retrieval of tubes



Pack of 4

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS23229B	Blue ■	2–5	50	12.9	7.8 x 4.1 x 0.9	19.9 x 10.5 x 2.8

**CORPORATE OFFICE
GENAXY SCIENTIFIC PVT. LTD**

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com