



Work 2 Store® Expanding Storage Box Single rack, double function

Streamline and reduce sample handling time with this convenient working and storage rack, eliminating the need several racks when working with samples requiring cold temperature storage

Efficient Handling

- Expand rack to opening position instantly creating more working space.
- Enough room for sample manipulation with no interference from open microfuge lids
- Collapse rack into closed position when it's ready for freezer storage
- Compact design will fit into standard 53 mm freezer racking system
- Alphanumeric identification for quick location and retrieval of tubes

Safe and Stackable Storage

- Clear lids locks the rack together for safe handling and transportation
- Low profile, easily stackable

Holds 64 x 1.5/2.0 mL cryo and/or microtubes

- Expanding wells have a center depression to stabilize conical tubes and a cryogenic vial lock for one handed cap removal



Work 2 Store® Expanding Tube Storage Box

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H (Open) in	L x W x H (Closed w/lid) in	L x W x H (Open) cm	L x W x H (Closed w/lid) cm
120018	White	1.5/2.0	64	12.3	12.2 x 5 x 1.3	5.2 x 5.2 x 2	31 x 12.8 x 3.4	13.3 x 13.3 x 5.3

CORPORATE OFFICE GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com