


TracRack® Tube Rack

Take what you need and store the rest

Comprised of a storage box and eight removable tube segments allowing you to easily take what you need - an entire segment or individual sample - and store the rest!

Convenient Rack System

- Eight polypropylene segments hold 1.5–2.0 mL cryo and microcentrifuge tubes (64 tubes)
- Wells feature a center depression stabilizing conical tubes and a cryogenic vial lock for one-handed cap removal
- Printed index letters on each segment make for quick sample ID and retrieval



TracRack® Outer Box

- Made from thin film polypropylene material
- Moisture and mold resistant
- 2 in (5.3 cm) boxes fit in standard freezer racks
- Not recommended for autoclaving

TracRack® Tube Rack

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120356	Purple	1.5/2.0	64	12.3	5.1 x 5.1 x 1.4	12.9 x 12.9 x 3.5

CORPORATE OFFICE
GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com