


Foam Tube Racks

High storage, minimal footprint

- Two-piece racks are constructed of expanded polystyrene (EPS) foam
- Lightweight microtube racks are ideal for long-term studies or tube storage
- Extra spacing between wells allows for easy gripping with fingertips
- Racks stack securely with nesting features on lid and base
- Must use tubes capable of being stored at the required temperature of the user



Pack of 100



Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS4282	White	1.5/2.0	100	11.5	13.2 x 3.7 x 2.9	33.6 x 9.5 x 7.3
HS4283P	White	1.5/2.0	50	11.5	8.3 x 4.3 x 2.7	21 x 11 x 7.1