

## Cryogenic Vial Boxes & Partitions

### For storage of small samples and tubes

- Made from cardboard
- Colors allow for easy identification and sorting purposes
- All box footprints are 13.3 x 13.3 cm, select from two different heights



### Cardboard Cryogenic Vial Boxes

Pack of 12


Item No.	Color		Tubes mL	L x W x H in	L x W x H cm
HS2860A	White	□	1.5/2.0	5.2 x 5.2 x 2	13.3 x 13.3 x 5
HS2860B	White	□	1.5/2.0	5.2 x 5.2 x 3	13.3 x 13.3 x 7.5
HS2860C	Assorted	■ ■ ■			
HS2860CB	Blue	■	1.5/2.0	5.2 x 5.2 x 2	13.3 x 13.3 x 5
HS2860CR	Red	■			
HS2860CY	Yellow	■			

### Partitions

Partitions must be ordered separately

Pack of 12

Item No.	Array	Tubes mL	Wells	Well Dia. mm
HS2860E	10 x 10	1.5/2.0	100	11.5
HS2860H	9 x 9	1.5/2.0	81	12
HS2860F	8 x 8	1.5/2.0	64	15
HS2860G	7 x 7	1.5/2.0	49	17.5



CORPORATE OFFICE  
GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058  
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com