


Microtube Storage Boxes

Safe microtube storage

Low profile boxes are available in two sizes to accommodate standard 0.5 mL and 0.2mL microtubes and vessels of a similar size

- Partitions included in the boxes
- Will fit into most standard freezer racking



Pack of 10

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120220	White <input type="checkbox"/>	0.2	169	9.5	5.2 x 5.2 x 1.1	13.3 x 13.3 x 2.9
120221	White <input type="checkbox"/>	0.5	100	12	5.2 x 5.2 x 1.6	13.3 x 13.3 x 4.1

CORPORATE OFFICE
GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com