

## Freezer Boxes Reusable and Affordable

Easy to assemble boxes made from polypropylene material

- Moisture resistant, more durable than cardboard boxes
- Dividers included with the boxes
- Corrugated Polypropylene: Store down to temperatures of -80°C with appropriate temperature viable tubes
- Thin Film: Can be used in mechanical freezers & vessels down to -196°C with appropriate temperature viable tubes
- Vapor phase only



Freezer boxes are supplied "flat-packed"  
Eliminating the hassle of stacking and reduces the  
amount of space needed for storage

### Corrugated Polypropylene 0.2/0.5/1.5–2.0/5 mL/Vial



Pack of 10

Item No.	Color		Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120379	White	□	0.2	144	7.9	5.2 x 5.1 x 1.2	13.3 x 13 x 3
120383	White	□	0.5	81	12	5.2 x 5.1 x 1.6	13.3 x 13 x 4.2
120262	White	□	1.5/2.0	81	12	5.2 x 5.1 x 2	13.3 x 13 x 5.1
120376	Blue	■	—	—	—	—	—
120387	White	□	5	25	22.6	5.2 x 5.1 x 2.9	13.3 x 13 x 7.5
120388	Blue	■	—	—	—	—	—
120428	White	□	Vial	100	17	7.4 x 7.8 x 1.6	18.8 x 19.8 x 4.1

### Corrugated Polypropylene 15 mL and 50 mL



Pack of 10

Item No.	Color		Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120391	White	□	15	36	19.8	5.8 x 5.7 x 4.7	14.7 x 14.5 x 12
120392	Blue	■	—	—	—	—	—
120366	White	□	50	16	31.6	5.8 x 5.7 x 4.7	14.7 x 14.5 x 12
120367	Blue	■	—	—	—	—	—

### Thin Polypropylene film 1.5–2.0 mL



Pack of 10

Item No.	Color		Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120370	White	□	—	—	—	—	—
120371	Blue	■	1.5/2.0	81	12.1	5 x 4.9 x 1.9	12.5 x 12.4 x 4.9
120372	Lt. Green	■	—	—	—	—	—

**CORPORATE OFFICE**  
**GENAXY SCIENTIFIC PVT. LTD**

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058  
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com