

Cargo Pipette Rack


Meets the need of any lab

Racks securely lock together horizontally or vertically to desired configuration or can be used individually

Can be used on top of benchtop or conveniently hung on any metallic surface with strong magnets that are designed into the racks

- Angled shelves store 0.1 to 50 mL pipettes
- Protective lid encloses the unit to prevent dust build-up and provides additional shelf space

Magnets
connect securely to
most metallic surfaces



Units Lock Securely
together horizontally or vertically

Pack of 3

Item No.	Color	L x W x H <i>in</i>	L x W x H <i>cm</i>
HS20615M	Blue ■	11.9 x 3.9 x 5	30.1 x 10 x 12.8

CORPORATE OFFICE
GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com