



Pipette Racks

Protects and organizes glass and plastic pipettes

Four angled-shelved compartments keep your 0.1 to 50 mL pipettes or cryo vial canes within easy reach while preventing damage and clutter

- Heavy-duty acrylic version provides crystal-clear viewing from any angle
- ABS versions have cut-out viewing window and come with two removable attachments that hold pipette bulbs, fillers, and pumps
- HS20612M magnetized ABS version has strong magnets on both sides and mounts securely to steel surfaces



Pipette Attachments

ABS versions includes
2 removable pipette attachments

Item No.	Color	Material	L x W x H <i>in</i>	L x W x H <i>cm</i>
HS20611	Clear <input type="checkbox"/>	Acrylic		
HS20612	Blue <input checked="" type="checkbox"/>	ABS	11.4 x 3.5 x 15.8	29.1 x 8.8 x 40.1
HS20613	Red <input checked="" type="checkbox"/>	ABS		
HS20612M	Blue <input checked="" type="checkbox"/>	ABS (mounting magnets)		

CORPORATE OFFICE
GENAXY SCIENTIFIC PVT. LTD

155-156, C-Block, Community Centre, Plot 32 to 57. Behind Janak Cinema, Janak Puri, New Delhi – 110058
Ph: +91-11-4761 9999(100 Lines) Web: www.genaxy.com E-Mail: technical@genaxy.com