



## 50 HOLE PC RACK

Material: Polycarbonate

50-Hole Rack can store up to fifty 1.2mL & 2.0mL tubes.

Chemical Resistant

Perfect for cryogenic storage

Color-coded for easier identification of samples

Stackable

The locking feature on the base enables speedy and convenient one-hand operation.

Cat. No.	Product Description	Pack Size
GEN-CVR-50-Y	50-Hole PC Rack for 1.2ml/2.0ml tubes, (L)143.6 x(W)73.8 x (H)52.9 mm	5

**Genaxy Scientific Pvt. Ltd .**