



PP CRYO BOX

Genaxy polypropylene cryo box are perfect for cryogenic storage of 1.8 ml or 2.0ml cryovials.

These racks are autoclavable and are ideal for use in benchtop area.

Options of 50 well, 81 well and 100 well are available.

Cat. No.	Product Description	Pack Size
GEN-CSR-50	50-Hole PP Cryobox, 80 x 62 x 54 mm	5
GEN-CSR-81	81-Hole PP Cryobox	5
GEN-CSR-100	100-Hole PP Cryobox 140 x 140 x 54mm	5

Genaxy Scientific Pvt. Ltd .