



## WORKSTATION RACK

Genaxy polypropylene workstation racks comes with unique locking features which makes the handling more convenient. Molded handles helps in scure handling and transport of samples. Saves space with their stackable feature.

Holds up to 50 External or Internal thread Cryogenic Vials

Locking feature enables speedy and convenient one-hand operation

Cat. No.	Product Description	Pack Size
GEN-CV-WR-50	50-Hole Workstation Rack for External or Internal Thread Cryogenic Vials, White, (W) 270mm x (D) 104mm x (H) 23mm	5

**Genaxy Scientific Pvt. Ltd .**