



96 HOLE PCR RACK WITH HINGED LID

Genaxy PCR racks with hinged lid can hold up to 96 samples which can be easily identified with alpha-numeric marking.

Samples are securely stored with lid

Holds up to ninety-six 0.2mL PCR tubes

Wells are numbered and lettered for easy identification

Cat. No.	Product Description	Pack Size
GEN-PCR-R96	96-Hole PCR Rack (with Hinged Lid) for 0.2ml PCR Tubes (L)127mm x (W)90mm x (H)31mm	5

Genaxy Scientific Pvt. Ltd .

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058
Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : info@genaxy.com