



96 HOLE TWO WAY MICROTUBE RACK

Genaxy two way microtube racks enables the use of two different volume in one rack. Cover lids prevents contamination.

Material : (Rack) Polypropylene / (Lid) Polycarbonate
Autoclavable

Used for 0.5mL, 1.5mL, and 2.0mL microtubes

Cover Lid prevents contamination

Rack Dimensions: (L)210mm x (W)110mm x (H)28mm

Cover Dimensions: (L)204mm x (W)115mm x (H)29mm

For use with 2.0mL Tube: Φ 11mm / Hole Depth 25.5mm / Flat bottom / 96 Hole

For use with 0.5mL Tube: Φ 7.8mm / Hole Depth 25.5mm / Flat bottom / 96 Hole

Cat. No.	Product Description	Pack Size
GEN-TWR-96	96-Hole Two-Way Microtube Rack for 0.5 ml & 1.5/2.0 ml Microtubes	5

Genaxy Scientific Pvt. Ltd .

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058
Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : info@genaxy.com