



## 80 HOLE TWO WAY MICROTUBE RACK

Genaxy two way microtube racks enables the use of two different volume in one rack. The racks are made up of polypropylene and are autoclavable. Comes without lid

Materials : Polypropylene

Autoclavable

Used for 0.5mL, 1.5mL, and 2.0mL microtubes

Dimensions: (L)225mm x (W)66mm x (H)26mm

For use with 2.0mL Tube:  $\Phi$ 10.8mm / Hole Depth 25mm / Flat bottom / 80 Hole

For use with 0.5mL Tube:  $\Phi$ 6.5mm / Hole Depth 19mm / Flat bottom / 60 Hole

| Cat. No.   | Product Description  | Pack Size |
|------------|--|-----------|
| GEN-TWR-80 | 80-Hole Two-Way Microtube Rack for 0.2ml, 1.5ml & 2.0ml Microtubes | 5         |

**Genaxy Scientific Pvt. Ltd .**

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : info@genaxy.com