



## CONICAL TUBE RACK

Genaxy Conical Tube Racks are manufactured with Polypropylene material and is autoclavable. These racks are very lightweighted and easy to assemble. Specific base design to hold conical tubes well at place. These can accommodate 12 numbers of 15ml and 15 numbers of 50ml tubes.

Material : Polypropylene

Accomodate:15mL, 50mL Conical Tubes

50mL Hole: 15 EA,  $\Phi$ 29.3

15mL Hole: 12 EA,  $\Phi$ 17

Dimensions: (L)233mm x (W)140mm x (H)48mm

Cat. No.	Product Description	Pack Size
GEN-PPR-1550	Racks for15mL and 50mL Conical Tubes	5

**Genaxy Scientific Pvt. Ltd .**

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058  
Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : info@genaxy.com