



## 1.5ML MICROTUBE RACK

Genaxy 1.5ml microtube racks are easy to assemble and lightweighted racks. These racks can be autoclaved.

Material : Polypropylene

32 place rack for 1.5mL Microtubes

Autoclavable

Lightweighted and easy to assemble

Dimensions : (L)158mm x (W)104mm x (H)30mm

Cat. No.	Product Description	Pack Size
GEN-PPR-20-32	2.0mL Microtube and Cryovial	5

**Genaxy Scientific Pvt. Ltd .**