



EPOXY COATED CENTRIFUGE TUBE

Genaxy epoxy coated racks are made of stainless steel. Epoxy coating prevents the racks from rust. These are autoclavable and easy to handle.

Material: Stainless Steel

Used for 15mL, 50mL and 250mL Centrifuge Tubes

Epoxy-coated to prevent rust

Autoclavable, and easy to handle

Dimensions (CWR-15) : (L)152mm x (W)82 mm x (H)83mm

Dimensions (CWR-50) : (L)151.5mm x (W)82.5mm x (H)82mm

Dimensions (CWR-100) : (L) 227mm x (W) 127mm x (H) 73mm

Dimensions (CWR-250) : (L)208mm x (W)136mm x (H)103mm

Dimensions (CWR-500) : (L) 210mm x (W) 205mm x (H) 100mm

Cat. No.

Product Description

Pack Size

GEN-CWR-15	18-Hole Epoxy Coated Centrifuge Tube Rack for 15ml Centrifuge Tubes, 18mm Hole Size	1
GEN-CWR-50	8-Hole Epoxy Coated Centrifuge Tube Rack for 50ml Centrifuge Tubes, 30mm Hole Size	1
GEN-CWR-100	8-Hole Epoxy Coated Centrifuge Tube Rack for 100ml Centrifuge Tubes, 52mm Hole Size	1
GEN-CWR-250	6-Hole Epoxy Coated Centrifuge Tube Rack for 250ml Centrifuge Tubes, 60mm Hole Size	1
GEN-CWR-500	4-Hole Epoxy Coated Centrifuge Tube Rack for 500ml Centrifuge Tubes, 97mm Hole Size	1

Genaxy Scientific Pvt. Ltd .

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : info@genaxy.com