



EVA FOAM RACK AND PP TUBE RACK

Genaxy Foam Racks is best suited for ice pans and buckets and can be used in water, melted ice and slurries.

1.5mL & 2.0mL PP Tube Rack holds up to 35 tubes

15mL PP Tube Rack holds up to 17 tubes

1.5mL & 2.0mL Square Foam Rack holds up to 63 tubes

15mL Square Foam Rack holds up to 35 tubes

Cat. No.	Product Description	Pack Size
GEN-IA-15FR	Round Foam Rack for 15mL tubes	1
GEN-IA-EFR	Round Foam Rack for 1.5mL & 2.0mL tubes	1
GEN-IA-15FS	35-Hole Square Foam Rack for 15mL tubes	1
GEN-IA-EFS	63-Hole Square Foam Rack for 1.5mL & 2.0mL tubes	1
GEN-IA-EP-W	35-Hole PP Tube Rack for 1.5mL & 2.0mL tubes	1
GEN-IA-15P-W	17-Hole PP Tube Rack for 15mL tubes	1

Genaxy Scientific Pvt. Ltd .