

Heavy Duty Vacuum Bottle, PP



- Polypropylene material
- Option for Medical Grade PP
- Medical Grade conforming to USP Class VI, US FDA CFR
- Meets Chemical Registration of TSCA (USA) & DSL (Canada).
- Autoclavable

Cat. Number	Description	Volume	Packaging	Type/Style	Qty. / Pack
GEN-HDPEVB-PP-2L	Heavy Duty Vacuum Bottle	2000ml	Bulk	PP	6
GEN-HDPEVB-PP-4L	Heavy Duty Vacuum Bottle	4000ml	Bulk	PP	6
GEN-HDPEVB-PP-2L-MG	Heavy Duty Vacuum Bottle	2000ml	Bulk	PP, Medical Grade	6
GEN-HDPEVB-PP-4L-MG	Heavy Duty Vacuum Bottle	4000ml	Bulk	PP, Medical Grade	6

Genaxy Scientific Pvt. Ltd.