

## Wide Mouth Bottle, Polypropylene



- Polypropylene Material
- Leak Proof
- Fully autoclavable

Cat. Number	Description	Volume	Packaging	Type/Style	Qty. / Pack
GEN-WMB-PP-30	Wide Mouth Bottle	30ml	Bulk	Polypropylene	72
GEN-WMB-PP-60	Wide Mouth Bottle	60ml	Bulk	Polypropylene	72
GEN-WMB-PP-125	Wide Mouth Bottle	125ml	Bulk	Polypropylene	72
GEN-WMB-PP-250	Wide Mouth Bottle	250ml	Bulk	Polypropylene	72
GEN-WMB-PP-500	Wide Mouth Bottle	500ml	Bulk	Polypropylene	48
GEN-WMB-PP-1L	Wide Mouth Bottle	1000ml	Bulk	Polypropylene	24

## Wide Mouth Bottle, LDPE



- Low Density Polyethylene (LDPE)
- Leak Proof
- Individually packed option available
- Non-autoclavable

Cat. Number	Description	Volume	Packaging	Type/Style	Qty. / Pack
GEN-WMB-LDPE-30	Wide Mouth Bottle	30ml	Bulk	LDPE	72
GEN-WMB-LDPE-60	Wide Mouth Bottle	60ml	Bulk	LDPE	72
GEN-WMB-LDPE-125	Wide Mouth Bottle	125ml	Bulk	LDPE	72
GEN-WMB-LDPE-250	Wide Mouth Bottle	250ml	Bulk	LDPE	72
GEN-WMB-LDPE-500	Wide Mouth Bottle	500ml	Bulk	LDPE	48
GEN-WMB-LDPE-1L	Wide Mouth Bottle	1 lt.	Bulk	LDPE	24
GEN-WMB-LDPE-30-S	Wide Mouth Bottle	30ml	Individually Pack	LDPE, Sterile	24
GEN-WMB-LDPE-60-S	Wide Mouth Bottle	60ml	Individually Pack	LDPE, Sterile	24
GEN-WMB-LDPE-125-S	Wide Mouth Bottle	125ml	Individually Pack	LDPE, Sterile	24
GEN-WMB-LDPE-250-S	Wide Mouth Bottl	250ml	Individually Pack	LDPE, Sterile	24
GEN-WMB-LDPE-500-S	Wide Mouth Bottle	500ml	Individually Pack	LDPE, Sterile	12
GEN-WMB-LDPE-1L-S	Wide Mouth Bottle	1 lt.	Individually Pack	LDPE, Sterile	12

# WIDE MOUTH BOTTLES

## Wide Mouth Bottle, HDPE



- High Density Polyethylene (HDPE) material
- Leak Proof
- Non-autoclavable

Cat. Number	Description	Volume	Packaging	Type/Style	Qty. / Pack
GEN-WMB-HDPE-30	Wide Mouth Bottle	30ml	Bulk	HDPE	72
GEN-WMB-HDPE-60	Wide Mouth Bottle	60ml	Bulk	HDPE	72
GEN-WMB-HDPE-125	Wide Mouth Bottle	125ml	Bulk	HDPE	72
GEN-WMB-HDPE-250	Wide Mouth Bottle	250ml	Bulk	HDPE	72
GEN-WMB-HDPE-500	Wide Mouth Bottle	500ml	Bulk	HDPE	48
GEN-WMB-HDPE-1L	Wide Mouth Bottle	1000ml	Bulk	HDPE	24

## Wide Mouth Bottle, HDPE, Amber



- High Density Polyethylene (HDPE) bottles
- Amber Color for light sensitive samples
- Opaque
- Non-autoclavable

Cat. Number	Description	Volume	Packaging	Type/Style	Qty. / Pack
GEN-WMB-HDPE-30-A	Wide Mouth Bottle	30ml	Bulk	HDPE Amber	72
GEN-WMB-HDPE-60-A	Wide Mouth Bottle	60ml	Bulk	HDPE Amber	72
GEN-WMB-HDPE-125-A	Wide Mouth Bottle	125ml	Bulk	HDPE Amber	72
GEN-WMB-HDPE-250-A	Wide Mouth Bottle	250ml	Bulk	HDPE Amber	72
GEN-WMB-HDPE-500-A	Wide Mouth Bottle	500ml	Bulk	HDPE Amber	48
GEN-WMB-HDPE-1L-A	Wide Mouth Bottle	1000ml	Bulk	HDPE Amber	24

# WIDE MOUTH BOTTLES

## Wide Mouth Bottle With Handle, Polypropylene



- Polypropylene material
- Option for Medical Grade PP
- Medical grade bottles are conforming to USP Class VI, US FDA CFR
- Meets Chemical Registration of TSCA (USA) & DSL (Canada).
- Autoclavable

Cat. Number	Description	Volume	Packaging Type/Style	Qty. / Pack
GEN-NMB-HDPE-4	Narrow Mouth bottle	4ml	Bulk	HDPE, AMBER 100
GEN-NMB-HDPE-8	Narrow Mouth bottle	8ml	Bulk	HDPE, AMBER 100
GEN-NMB-HDPE-15	Narrow Mouth bottle	15ml	Bulk	HDPE, AMBER 100
GEN-NMB-HDPE-30	Narrow Mouth bottle	30ml	Bulk	HDPE, AMBER 72
GEN-NMB-HDPE-60	Narrow Mouth bottle	60ml	Bulk	HDPE, AMBER 72
GEN-NMB-HDPE-125	Narrow Mouth bottle	125ml	Bulk	HDPE, AMBER 72
GEN-NMB-HDPE-250	Narrow Mouth bottle	250ml	Bulk	HDPE, AMBER 72
GEN-NMB-HDPE-500	Narrow Mouth bottle	500ml	Bulk	HDPE, AMBER 48
GEN-NMB-HDPE-1L	Narrow Mouth bottle	1000ml	Bulk	HDPE, AMBER 24

## Genaxy Scientific Pvt. Ltd.