

Wash Bottle



- Moulded from high quality virgin low-density polyethylene (LDPE) to ensure excellent results
- Smooth dispensing with ergonomic soft squeeze
- These wash bottles are used for routine storage and pouring applications in the laboratory
- Available in Bulk or Individually wrapped packing

Cat. Number	Description	Volume	Type/Style	Qty. / Pack
GEN-WB-LDPE-250	Wash Bottle	250ml	LDPE	12
GEN-WB-LDPE-500	Wash Bottle	500ml	LDPE	6
GEN-WB-LDPE-1L	Wash Bottle	1000ml	LDPE	6
GEN-WBIP-LDPE-125	Wash Bottle	125ml	Individually Wrapped, LDPE	12
GEN-WBIP-LDPE-250	Wash Bottle	250ml	Individually Wrapped, LDPE	12
GEN-WBIP-LDPE-500	Wash Bottle	500ml	Individually Wrapped, LDPE	6

Genaxy Scientific Pvt. Ltd.