

Measuring Jug, PP



- The Measuring Jugs, moulded in Polypropylene are clear
- Autoclavable
- Good chemical resistance
- Non-drip spout
- Provided with a handle with excellent thumb grip for easy & convenient lifting
- Raised graduations for easy reading

| Cat. Number | Description | Volume | Type/Style | Qty. / Pack |
|---------------|---------------|--------|------------|-------------|
| GEN-MJ-PP-250 | Measuring Jug | 250ml | PP | 12 |
| GEN-MJ-PP-500 | Measuring Jug | 500ml | PP | 6 |
| GEN-MJ-PP-1L | Measuring Jug | 1000ml | PP | 6 |
| GEN-MJ-PP-2L | Measuring Jug | 2000ml | PP | 2 |

Genaxy Scientific Pvt. Ltd.