

Beakers, Euro, PP



- Beakers, moulded in Polypropylene
- Excellent clarity
- Very good chemical resistance
- Autoclavable
- Prominently raised graduations on the outer surface
- Pouring a Solution becomes very convenient because of its tapered spout

Cat. Number	Description	Volume	Type/Style	Qty. / Pack
GEN-BK-EPP-50	Beaker	50ml	Euro - PP	12
GEN-BK-EPP-100	Beaker	100ml	Euro - PP	12
GEN-BK-EPP-250	Beaker	250ml	Euro - PP	12
GEN-BK-EPP-500	Beaker	500ml	Euro - PP	12
GEN-BK-EPP-1L	Beaker	1000ml	Euro - PP	6
GEN-BK-EPP-2L	Beaker	2000ml	Euro - PP	2
GEN-BK-EPP-5L	Beaker	5000ml	Euro - PP	1

Genaxy Scientific Pvt. Ltd.