

Funnel with Long Stem, PP



- Polypropylene (PP) analytical funnels are perfect for use in analytical Chemistry applications
- Features long stems making them ideal for use with liquids
- Internal and external ribbing to prevent air lock.
- Angled at exactly 60° for rapid
- Filtration for use with standard filter paper
- Autoclavable at 121°C for 20 minutes

| Cat. Number | Description | Volume | Type/Style | Qty. / Pack |
|----------------|-----------------------|--------|------------|-------------|
| GEN-FLS-PP-35 | Funnel for Long, Stem | 35ml | PP | 36 |
| GEN-FLS-PP-50 | Funnel for Long, Stem | 50ml | PP | 36 |
| GEN-FLS-PP-65 | Funnel for Long, Stem | 65ml | PP | 36 |
| GEN-FLS-PP-75 | Funnel for Long, Stem | 75ml | PP | 36 |
| GEN-FLS-PP-100 | Funnel for Long, Stem | 100ml | PP | 12 |

Genaxy Scientific Pvt. Ltd.