

The perfect ECL formulation combining great sensitivity and long signal duration



WESTAR ANTARES is an ECL substrate with stable light output and mid-femtogram detection level

Cat#	Description	Pack Size	Qty	Unit Price	Offer*
XLS0142,0250	WESTAR ANTARES	2 x 125 mL kit	1 Kit	16,800/-	
	+				
XLS0142,0250	WESTAR ANTARES	2 x 125 mL kit	1 Kit	16,800/-	16,800/-

Features

Mid femtogram detection
Extended signal duration
Working solution stable for at least three days
High range flexibility

Components

Luminol /enhancer solution (A) 125ml and Peroxide solution (B) 125ml (for 2500 cm² of membrane).

Antibody Dilution Range

Primary: 1:1.000 - 1:15.000
Secondary: 1:25.000 - 1:150.000 (from 1 mg/mL stock solution)

Signal duration

> 24 hours

Storage

Stable for 1 year at RT



***Terms and Condition**

Valid till 31st July 2022

Govt. Tax as per applicable will be charged extra

This offer can not be clubbed with any other offer or RC discount

GENAXY SCIENTIFIC PRIVATE LIMITED

155-156, C Block, Community centre Plot No. 32 to 57, Behind Janak Cinema Janakpuri, New Delhi - 110058, Email: technical@genaxy.com, Ph : +91-11-47619999 (100 Lines), Mob:9971027408