



4-WAY UNIVERSAL TUBE RACKS

Heavy-duty polypropylene racks

Versatile design accommodates different tube sizes on each side

The four sides holds 50ml, 15ml, 1.5/2.0ml and 0.5ml tubes. Racks are autoclavable and available in different colors.

Cat. No.	Product Description	Pack Size
GEN-UTR	4-way Universal Tube Rack for 4 x 50ml tubes / 12 x 15ml tubes / 32 x 0.5ml tubes / 32 x 1.5 or 2.0ml tubes, PP, Autoclavable, 17.4 x 9.5 x 5.2 mm	4

Genaxy Scientific Pvt. Ltd .

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058
Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : info@genaxy.com