

GENAXY Mini Centrifuge

GEN-MINI-6KC

Mini-6KC Mini-Centrifuge is compact and cost effective, allowing any workstation to be equipped with a “personal” mini centrifuge. It includes 2 quick release interchangeable rotors, two speed options (4,000rpm & 6000rpm). The electronic brake can reach quick deceleration to reduce handling time.

Technical Specification

Model	GEN-MINI-6KC
Speed	4,000rpm & 6,000rpm
Max. Capacity	8 x 1.5/2.0 ml tube rotor 4 place PCR Strip rotor
Number of Rotors	2
Maximum RCF	1,200 & 2,000 x g
Dimensions(D×W×H)	135 x 147 x 107mm
Weight	0.7kg
Power Supply	AC100~240V, 50/60Hz



GEN-MINI-6KC



Features

Composite shell, shatter-resistant

2 quick release rotors are available

Electronic braking system for rapid processing

Max. speed: 6000 rpm (2000 × g)

Quick release rotor tub

easy rotor exchange

8x1.5/2.0ml tube rotor, 4 place PCR Strip rotor,

6x1.5/2.0ml tubes, --6x0.5ml tubes, --2-8x0.2ml PCR strip tubese

Accessories

S No.	Code	Description
1	AS-08021-01	8x1.5/2.0ml rotor
2	AS-08021-02	4x8x0.2ml PCR strip rotor
3	AS-08061-05	Adapters for 0.5/0.6ml tubes
4	AS-08061-06	Adapters for 0.2ml tubes

Genaxy Scientific Pvt. Ltd :

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163

E-mail : info@genaxy.com